**INTENDED USE**

Remel Xpect® Clostridium difficile Toxin A/B test kit is a rapid *in vitro* immunochromatographic test for the direct, qualitative detection of *Clostridium difficile* Toxin A and/or B in human fecal specimens from patients suspected of having *Clostridium difficile*-associated disease (CDAD). The test is intended for use as an aid in diagnosis of CDAD. The test can also be used for confirmation of toxigenic *C. difficile* from Brain Heart Infusion (BHI) broth culture.

SUMMARY AND EXPLANATION

C. difficile was first described in 1935, but it was not until 1977 that links were made between the organism and disease. *C. difficile* is a toxin-producing, spore-forming anaerobic gram-positive bacillus. The clinical presentation of *C. difficile* infection includes, in increasing order of severity, asymptomatic carriers, antibiotic-associated colitis, pseudomembranous colitis (PMC), and fulminant colitis.¹ A reduction of the normal microbial flora in the colon, usually caused by antibiotic therapy, allows overgrowth of *C. difficile*. Symptoms of antibiotic-associated colitis usually begin four to ten days after antibiotic treatment has begun.

Most pathogenic strains of *C. difficile* produce two toxins, Toxin A (enterotoxin) and Toxin B (cytotoxin), which are the main virulence factors for the organism. There appears to be a cascade of events, which result in the expression of the activity of these toxins. Toxin A is mildly cytopathic but induces large fluid shifts and mucosal inflammation. Toxin B is intensely cytopathic but its role in the disease process is not clearly understood. Variant strains which are Toxin A-negative, Toxin B-positive are known to exist, are fully pathogenic, and capable of producing the full spectrum of disease. The prevalence of these variant strains varies widely by institution and geographic location.²⁻⁷

C. difficile-associated disease (CDAD) primarily occurs in hospitalized patients.^{8,9} Individuals with CDAD shed spores in the stool, which can survive for as long as five months in the environment. Clinical and pathological features of CDAD are not easily distinguished from those of other gastrointestinal diseases, including ulcerative colitis, chronic inflammatory bowel disease, and Crohn's disease. Infection with toxigenic *C. difficile* is a potentially life-threatening disease process; however, when properly treated, patient mortality rates are low. Thus, rapid diagnosis, allowing clinicians to initiate appropriate therapy and implement adequate measures to control nosocomial spread, is important.

PRINCIPLE

The Xpect® Clostridium difficile Toxin A/B test is a qualitative immunochromatographic assay that detects *C. difficile* Toxin A and Toxin B in stool specimens or cultures of toxigenic *C. difficile*. In performing the test, a specimen is first diluted with Specimen Diluent to help solubilize the toxins. A portion of the diluted sample is then mixed with a volume of Conjugate 1 containing antibodies to Toxin A and Toxin B coupled to colored microparticles, plus a volume of Conjugate 2 containing biotinylated antibodies to Toxin A and Toxin B. A volume of this mixture is transferred to a test device having immobilized streptavidin as a test line and goat anti-immunoglobulin antibody as a control line. Immunocomplexes of toxin and conjugated antibodies form a visible band as they flow across the test line. Excess colored particle conjugates form a visible band at the control line to document that the test is functioning properly.

STORAGE

Store test devices in sealed foil pouches at 2-30°C (room temperature or refrigerated). Store all kit reagent bottles and vials at 2-8°C. Do not freeze or overheat. Allow components to equilibrate to room temperature before use. Mix bottled reagents gently prior to use. Return the unused reagents to the refrigerator after use.

PRECAUTIONS

1. For *In Vitro* Diagnostic Use.
2. Standard precautions should be taken against the dangers of biological hazards by properly sterilizing specimens, containers, and test devices after use. Consult appropriate references when necessary.¹⁰
3. Directions should be read and followed carefully.
4. Reagents are provided at the necessary working strength and are to be dispensed directly from the dropper bottles. Do not dilute reagents.
5. Do not interchange reagents between kits of different lots.
6. Do not use reagents beyond the printed expiration dates.
7. Microbial contamination of reagents may decrease the accuracy of the assay.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected in clean, airtight, leak-proof containers.

- Fresh, untreated stool specimens should be stored at 2-8°C and tested within 72 hours of collection. If fresh specimens cannot be tested within 72 hours, they should be frozen at -20°C or below in a non-defrosting freezer and tested within 2 months of collection. Avoid multiple freeze-thaw cycles.
- Fresh specimens diluted in the Sample Diluent provided in the kit can be stored refrigerated for up to 24 hours prior to testing.
- Stool specimens collected in modified Cary Blair Transport Medium with indicator (or equivalent) may be stored refrigerated (2-8°C) or stored at room temperature (20-25°C) and should be tested within 5 days of collection.
- Stool specimens that have been concentrated or collected in Formalin, SAF, or PVA are not suitable for use with this test.

REAGENTS AND MATERIALS SUPPLIED**20 / 40 Tests**

1. **Test Devices** (20 / 40): Each foil pouch contains one single-use test device with desiccant; membrane is striped with capture reagents
2. **Conjugate Reagent 1** (3.0 ml x 1) / (3.0 ml x 2): Blue-black microparticles coated with mouse anti-toxin A and rabbit anti-Toxin B
3. **Conjugate Reagent 2** (3.0 ml x 1) / (3.0 ml x 2): Biotinylated goat anti-Toxin A and rabbit anti-Toxin B
4. **Specimen Diluent** (12.0 ml x 1) / (12.0 ml x 2): Buffered solution with 0.03% ProClin™ 300
5. **Positive Control** (2.0 ml x 1) / (2.0 ml x 2): Culture supernatant containing *C. difficile* Toxin A and B with preservative
6. **Negative Control** (2.0 ml x 1) / (2.0 ml x 2): Buffered solution with preservative
7. **Disposable Transfer Pipettes** (40 / 80): Pipettes with marked graduations at 0.1 ml increments
8. **Dilution Tubes** (40 / 80): Tubes for specimen preparation with marked 0.5 ml volume
9. **Wooden Applicator Sticks** (20 / 40)
10. **Instructions for Use** (1) (IFU)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Stool specimen collection container(s), and (2) Timer.

Optional materials not provided:

(1) Specimen transport media, (2) BHI Broth, and (3) Test tube rack.

PROCEDURE

1. Allow kit components and stool specimens to come to room temperature before use.
2. Mix stool specimens **thoroughly** prior to testing (regardless of consistency).
3. Remove the test device from the foil pouch when ready to perform the test and place it on a flat surface.
4. Label the device (and dilution tubes) with patient or control identification.

Sample Preparation

Fresh (unpreserved) Stool Specimens:

1. Add Specimen Diluent up to the line marked on the dilution tube (0.5 ml).
 - a. For solid stool specimens, use a wooden applicator stick and sample from separate areas of the specimen. Transfer a 0.2 g (6-8 mm diameter) portion of stool to the dilution tube.
 - b. For semi-solid or liquid stool specimens, use a transfer pipette to transfer 0.2 ml (second graduated mark from the tip) of specimen to the dilution tube.
2. Mix thoroughly and allow large particulates to settle.
3. Use a transfer pipette to dispense 0.1 ml (first molded graduated mark from the tip) of this preparation into another dilution tube.

Stool Specimens in Cary Blair:

1. Specimens collected in modified Cary Blair transport media with indicator are used without further dilution.
2. Mix specimen thoroughly prior to testing.
3. Use a transfer pipette to dispense 0.1 ml (first molded graduated mark from the tip) of the specimen into a dilution tube.

BHI Broth Culture:

1. Use a transfer pipette to dispense 0.1 ml (first molded graduated mark from the tip) of a 72-hour BHI broth culture with suspected *C. difficile* into a dilution tube.


Test Procedure

1. Add 5 drops of Conjugate Reagent 1 to the dilution tube containing the diluted specimen.
2. Add 5 drops of Conjugate Reagent 2 to the dilution tube containing the diluted specimen.
3. Mix the contents of the tube thoroughly.
4. Use a transfer pipette to dispense 0.2 ml (second graduated mark from the tip) of sample into the circular sample well of the test device.
5. Read and record the test results visually after 20 minutes according to the INTERPRETATION section. (Strong positive results may be apparent sooner than 20 minutes.) Results are invalid beyond 20 minutes.

INTERPRETATION OF THE TEST


Positive Result (Toxin A and/or B present):

A positive test is indicated by two black-colored lines of any intensity; one in the TEST region and one in the control (CTRL) region. A positive test indicates the presence of *C. difficile* Toxin A and/or Toxin B in the sample.


Negative Result (Toxin A and/or B not detected):

A negative test is indicated by only one black-colored line in the control (CTRL) region. A negative test indicates that *C. difficile* Toxin A and/or Toxin B are absent or below the detection limit of the test.

Invalid Result:

An invalid result occurs when the TEST line is partial or incomplete, or the control (CTRL) line is absent or incomplete at 20 minutes.

If invalid results occur due to restricted flow of the sample along the membrane, the specimen should be re-tested as follows:

Re-test Specimen Preparation

Further dilute the specimen as follows:

1. Combine in a clean dilution tube:
 - a. 4 drops of Specimen Diluent, and
 - b. 0.1 ml (first molded graduated mark from the tip) of specimen from the initial specimen preparation in Specimen Diluent (without reagents added).
2. Mix the preparation by aspirating with the transfer pipette.
3. Use 0.1 ml of this diluted specimen preparation and repeat the Test Procedure.

QUALITY CONTROL

Internal: A procedural control is included in the test. The appearance of a control line in the CTRL region verifies that an intact conjugate has been added to the device, that the control line antibody is functionally active, and that adequate capillary flow has occurred. A clear background in the results area is considered an internal negative control. If the test has been performed correctly and reagents are working properly, the background will clear to give a discernible result.

External: The Positive and Negative Quality Controls provided with the kit should be run with each shipment and new kit lot number received. The Positive Control is used to verify reactivity of the reagents associated with the assay and is not intended to ensure precision at the analytical assay cut-off. Each laboratory should follow their state and local requirements. To use, add 0.1 ml (5 drops) of the control to a dilution tube and process in accordance with Test Procedure outlined above (do not dilute controls with specimen diluent). If aberrant quality control results are noted, patient results should not be reported.

LIMITATIONS

1. A positive test does not define the presence of disease. The test detects the presence of Toxin A and/or Toxin B in fecal specimens. Results should be used in conjunction with other clinical findings to establish a diagnosis.
2. As with all *in vitro* diagnostic tests, a negative test result does not exclude the possibility of the presence of *C. difficile* Toxin A or Toxin B. This may occur when the toxin level in the sample is below the detection level of the test. The level of toxin has not been shown to correlate with either the presence or severity of disease. Test results should be interpreted by a physician in light of other laboratory results and clinical findings.
3. *Clostridium sordellii* is not a common inhabitant of the human intestine but does produce a toxin (HT) that has biological, physiochemical, and immunochemical properties similar to *C. difficile* Toxin A. This similarity may cause cross-reactivity in any diagnostic test that detects *C. difficile* Toxin A and/or B.¹¹

- Proper specimen collection and processing are essential to achieve optimal performance of the test. See Specimen Collection, Storage and Transport section.
- The performance characteristics of this test have not been fully evaluated in pediatric populations.
- A high dose Hook Effect was not observed during the BHI broth culture study in which *C. difficile* was evaluated after growing in pure culture for 72-hours.

EXPECTED VALUES

Community-acquired CDAD cases are recognized, but the incidence is low (<1 case per 10,000 antibiotic prescriptions). This may be because diagnostic testing is not performed often enough in outpatient settings to detect CDAD. *C. difficile* colitis occurs at a much higher frequency in patients who are hospitalized and is the fourth most common nosocomial disease reported to the Centers for Disease Control and Prevention. *C. difficile* is responsible for 20-30% of antibiotic-associated diarrhea and more than 90% of pseudomembranous colitis. The incidence rate of nosocomial CDAD may vary with hospital populations and is influenced by the presence of predisposing factors, such as increased patient age, type and duration of antimicrobial therapy, severity of underlying illness(es), and length of hospital stay. *C. difficile* is found in 3-5% of healthy adults and up to 50% of infants and young adults that asymptotically carry both the bacteria and its toxins.¹² An overall prevalence rate of 16% was observed when the Xpect® Clostridium difficile Toxin A/B was evaluated in a prospective study conducted at four independent laboratories in the United States.

PERFORMANCE CHARACTERISTICS

Clinical Accuracy:

The performance of Xpect® Clostridium difficile Toxin A/B was evaluated at four geographically diverse regions of the United States. A total of eight hundred fifteen specimens were tested with the Xpect® Clostridium difficile Toxin A/B test and compared to results obtained from the cytotoxin assay (CTA).

		CTA Results	
		+	-
Overall Xpect® Results	+	132	25
	-	21	637
TOTAL		153	662

Sensitivity: 86.3% (95% CI = 79.8-91.3%)
 Specificity: 96.2% (95% CI = 94.5-97.5%)
 Positive Predictive Value: 84.1% (95% CI = 77.4-89.4%)
 Negative Predictive Value: 96.8% (95% CI = 95.2-98.0%)
 % Correlation: 94.4% (95% CI = 92.5-95.8%)

Note: CI = Confidence Interval

Discordant results were further investigated by toxigenic culture and microwell enzyme immunoassay that detects both Toxin A and B. Four of 25 specimens that were cytotoxin negative and Xpect® Clostridium difficile Toxin A/B positive on initial testing were positive by toxigenic culture and enzyme immunoassay. Ten of 21 (47.6%) specimens that were CTA positive and Xpect® Clostridium difficile Toxin A/B negative on initial testing were negative by toxigenic culture and microwell enzyme immunoassay.

Performance Compared to Commercially Available Devices:

The Xpect® Clostridium difficile Toxin A/B test was also compared to two commercially available products. Each of the four clinical trial sites tested a chromatographic membrane assay that detects Toxin A only (Predicate Device), the Xpect® Clostridium difficile Toxin A/B test, and cytotoxin assay (CTA) for each sample. In addition, one clinical trial site tested a microwell enzyme immunoassay for the detection of both Toxin A and B for each sample. The results presented below are calculated using CTA as the reference.

Performance of Devices Compared to CTA:

	n = 815		n = 267	
	Xpect® C. difficile Toxin A/B	Predicate Device	Xpect® C. difficile Toxin A/B	EIA
Sensitivity	86.3%	62.7%	91.0%	80.6%
Specificity	96.2%	98.8%	98.0%	97.5%

BHI Broth Culture Performance:

An in-house study was conducted using twenty-one known reference strains and thirty-six suspect *C. difficile* isolates from stool specimens. BHI broth cultures were tested with Xpect® Clostridium difficile Toxin A/B following 72-hours incubation. Under these conditions, the BHI broth culture of *Clostridium sordellii* ATCC® 9714 produced a positive reaction. There was 94.7% (54/57) agreement with expected values.

Analytical Sensitivity:

The analytical sensitivity was evaluated using purified *C. difficile* Toxin A and Toxin B. The Xpect® Clostridium difficile Toxin A/B test detects Toxin A at levels of ≥ 6.25 ng/ml (0.12 ng/test) and Toxin B at levels of ≥ 40.0 ng/ml (0.76 ng/test).

Cross-Reactivity:

Fifty-four microorganisms were evaluated with Xpect® Clostridium difficile Toxin A/B. No cross-reactivity was observed. Bacteria and yeast isolates were tested at 10⁸ colony-forming units per ml concentration. Viral isolates were tested at concentrations of 10⁴ to 10⁵ TCID₅₀ (tissue culture infectious dose) per ml concentration. The following organisms were tested in the Xpect® Clostridium difficile Toxin A/B test:

<i>Aeromonas hydrophila</i>	<i>Klebsiella pneumoniae</i>
<i>Bacillus cereus</i>	<i>Peptostreptococcus anaerobius</i>
<i>Bacillus subtilis</i>	<i>Porphyromonas asaccharolytica</i>
<i>Bacteroides fragilis</i>	<i>Proteus mirabilis</i>
<i>Campylobacter coli</i>	<i>Proteus vulgaris</i>
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	<i>Salmonella</i> Typhimurium
<i>Campylobacter lari</i>	<i>Serratia liquefaciens</i>
<i>Candida albicans</i>	<i>Shigella boydii</i>
<i>Clostridium botulinum</i> (toxin 20 µg/ml)	<i>Shigella dysenteriae</i>
<i>Clostridium beijerinckii</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i> (non-toxigenic)	<i>Shigella sonnei</i>
<i>Clostridium haemolyticum</i>	<i>Staphylococcus aureus</i> (Cowan)
<i>Clostridium histolyticum</i>	<i>Staphylococcus epidermidis</i>
<i>Clostridium innocuum</i>	<i>Vibrio cholerae</i>
<i>Clostridium novyi</i>	<i>Vibrio parahaemolyticus</i>
<i>Clostridium perfringens</i>	<i>Yersinia enterocolitica</i>
<i>Clostridium septicum</i>	<i>Giardia intestinalis</i>
<i>Clostridium sordellii</i>	<i>Entamoeba histolytica</i>
<i>Clostridium sporogenes</i>	Adenovirus type 2
<i>Clostridium subterminale</i>	Adenovirus type 40
<i>Clostridium tetani</i>	Adenovirus type 41
<i>Enterobacter aerogenes</i>	Coxsackievirus B4
<i>Enterobacter cloacae</i>	Cytomegalovirus
<i>Enterococcus faecalis</i>	Echovirus (type 22)
<i>Enterococcus faecium</i>	Enterovirus (type 69)
<i>Escherichia coli</i>	Rotavirus

Interfering Substances:

The following substances were tested with the Xpect® Clostridium difficile Toxin A/B test and no interference was observed in the assay for any substance tested at the indicated levels: blood, mucous, fecal fat, Pepto-Bismol® (10%v/v), Imodium® AD (10%v/v), Kaopectate® (10%v/v), Castoria® (10%v/v), vancomycin (12.5 mg/ml), metronidazole (12.5 mg/ml), and barium sulfate (12.5 mg/ml).

Reproducibility:

Reproducibility testing was conducted at three sites, including one in-house site, on four separate days with six blinded samples. The samples consisted of known positive and negative stool specimens. The samples produced the expected results with the Xpect® Clostridium difficile Toxin A/B test 98.6% (71/72) of the time.

BIBLIOGRAPHY

1. Gerding, D.N., S. Johnson, L.R. Peterson, M.E. Mulligan, and J. Silva, Jr. 1995. Infect. Control Hosp. Epidemiol. 16:459-477.
2. Limaye, A.P., D.K. Turgeon, B.T. Cookson, and T.R. Fritzsche. 2000. J. Clin. Microbiol. 38:1696-1697.
3. Alfa, M.J., A. Kabani, D. Lyerly, S. Moncrief, L.M. Neville, A. Al-Barrak, G.K.H. Harding, B. Dyck, K. Olekson, and J.M. Embil. 2000. J. Clin. Microbiol. 38:2706-2714.
4. Barbut, F., V. Lalande, B. Burghoffer, H.V. Thien, E. Grimprel, and J. Petit. 2002. J. Clin. Microbiol. 40:2079-2083.
5. Johnson, S., S.P. Sambol, J.S. Brazier, M. Delmée, V. Avesani, M.M. Merrigan, and D.N. Gerding. 2003. J. Clin. Microbiol. 41:1543-1547.
6. Kato, H., N. Kato, K. Watanabe, N. Iwai, H. Nakamura, T. Yamamoto, K. Suzuki, S. Kim, Y. Chong, and E.B. Wasito. 1998. J. Clin. Microbiol. 36:2178-2182.
7. Johnson, S., S.A. Kent, K.J. O'Leary, M.M. Merrigan, S.P. Sambol, L.R. Peterson, and D.N. Gerding. 2001. Ann. Intern. Med. 135:434-438.
8. Wilkins, T.D. and D.M. Lyerly. 2003. J. Clin. Microbiol. 41:531-534.
9. Fedorka, D.P. 2002. Clin. Microbiol. Newsletter. 24:76-79.
10. Richmond, J.Y. and R.W. McKinney. 1999. Biosafety in Microbiological and Biomedical Laboratories. 4th ed. U.S. Government Printing Office. Washington, D.C.
11. Lyerly, D.M., H.C. Krivan, and T.D. Wilkins. 1988. Clin. Microbiol. Rev. 1:1-18.
12. Mandell, G.L., J.E. Bennett, and R. Dolin. 2000. Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases. 5th ed. Churchill Livingstone. New York, NY.

PACKAGING

REF R24650, Xpect® Clostridium difficile Toxin A/B 20 Tests/Kit
 REF R24640, Xpect® Clostridium difficile Toxin A/B 40 Tests/Kit

Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

Xpect® is a registered trademark of Remel Inc.

ATCC® is a registered trademark of American Type Culture Collection.
 Pepto-Bismol®, Imodium®, Kaopectate®, and Castoria® are registered trademarks of Procter & Gamble, McNeil Consumer & Specialty Pharmaceuticals, Pharmacia & Upjohn Company, and The Mentholatum Company, Inc., respectively.
 ProClin™ is a trademark of the Rohm and Haas Company.

IFU 24650, Revised August 27, 2007


Printed in U.S.A.

12076 Santa Fe Drive, Lenexa, KS 66215, USA

General Information: (800) 255-6730 Technical Service: (800) 447-3641 Order Entry: (800) 447-3635

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- Gerding, D.N., S. Johnson, L.R. Peterson, M.E. Mulligan, and J. Silva, Jr. 1995. Infect. Control Hosp. Epidemiol. 16:459-477.
- Limaye, A.P., D.K. Turgeon, B.T. Cookson, and T.R. Fritsche. 2000. J. Clin. Microbiol. 38:1696-1697.
- Alfa, M.J., A. Kabani, D. Lierly, S. Moncrief, L.M. Neville, A. Al-Barrak, G.K.H. Harding, B. Dyck, K. Olekson, and J.M. Embil. 2000. J. Clin. Microbiol. 38:2706-2714.
- Barbut, F., V. Lalande, B. Burghoffer, H.V. Thien, E. Grimprel, and J. Petit. 2002. J. Clin. Microbiol. 40:2079-2083.
- Johnson, S., S.P. Sambol, J.S. Brazier, M. Delmée, V. Avesani, M.M. Merrigan, and D.N. Gerding. 2003. J. Clin. Microbiol. 41:1543-1547.
- Kato, H., N. Kato, K. Watanabe, N. Iwai, H. Nakamura, T. Yamamoto, K. Suzuki, S. Kim, Y. Chong, and E.B. Wasito. 1998. J. Clin. Microbiol. 36:2178-2182.
- Johnson, S., S.A. Kent, K.J. O'Leary, M.M. Merrigan, S.P. Sambol, L.R. Peterson, and D.N. Gerding. 2001. Ann. Intern. Med. 135:434-438.
- Wilkins, T.D. and D.M. Lierly. 2003. J. Clin. Microbiol. 41:531-534.
- Fedorka, D.P. 2002. Clin. Microbiol. Newsletter. 24:76-79.
- Richmond, J.Y. and R.W. McKinney. 1999. Biosafety in Microbiological and Biomedical Laboratories. 4th ed. U.S. Government Printing Office. Washington, D.C.
- Lierly, D.M., H.C. Krivan, and T.D. Wilkins. 1988. Clin. Microbiol. Rev. 1:1-18.
- Mandell, G.L., J.E. Bennett, and R. Dolin. 2000. Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases. 5th ed. Churchill Livingstone. New York, NY.

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REF R24650, Xpect® Clostridium difficile Toxin A/B20 Tests/Kit
REF R24640, Xpect® Clostridium difficile Toxin A/B40 Tests/Kit

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- Barbut, F., V. Lalande, B. Burghoffer, H.V. Thien, E. Grimpel, and J. Petit. 2002. J. Clin. Microbiol. 40:2079-2083.
- Johnson, S., S.P. Sambol, J.S. Brazier, M. Delmée, V. Avesani, M.M. Merrigan, and D.N. Gerding. 2003. J. Clin. Microbiol. 41:1543-1547.
- Kato, H., N. Kato, K. Watanabe, N. Iwai, H. Nakamura, T. Yamamoto, K. Suzuki, S. Kim, Y. Chong, and E.B. Wasito. 1998. J. Clin. Microbiol. 36:2178-2182.
- Johnson, S., S.A. Kent, K.J. O'Leary, M.M. Merrigan, S.P. Sambol, L.R. Peterson, and D.N. Gerding. 2001. Ann. Intern. Med. 135:434-438.
- Wilkins, T.D. and D.M. Lyster. 2003. J. Clin. Microbiol. 41:531-534.
- Fedorka, D.P. 2002. Clin. Microbiol. Newsletter. 24:76-79.
- Richmond, J.Y. and R.W. McKinney. 1999. Biosafety in Microbiological and Biomedical Laboratories. 4th ed. U.S. Government Printing Office. Washington, D.C.
- Lyster, D.M., H.C. Krivan, and T.D. Wilkins. 1988. Clin. Microbiol. Rev. 1:1-18.
- Mandell, G.L., J.E. Bennett, and R. Dolin. 2000. Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases. 5th ed. Churchill Livingstone. New York, NY.

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REF R24650, Xpect® Toxine A/B de Clostridium difficile..... 20 tests/kit
REF R24640, Xpect® Toxine A/B de Clostridium difficile..... 40 tests/kit

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

BIBLIOGRAFIA

- Gerding, D.N., S. Johnson, L.R. Peterson, M.E. Mulligan, and J. Silva, Jr. 1995. Infect. Control Hosp. Epidemiol. 16:459-477.
- Limaye, A.P., D.K. Turgeon, B.T. Cookson, and T.R. Fritsche. 2000. J. Clin. Microbiol. 38:1696-1697.
- Alfa, M.J., A. Kabani, D. Lyerly, S. Moncrief, L.M. Neville, A. Al-Barrak, G.K.H. Harding, B. Dyck, K. Olekson, and J.M. Embil. 2000. J. Clin. Microbiol. 38:2706-2714.
- Barbut, F., V. Lalande, B. Burghoffer, H.V. Thien, E. Grimprel, and J. Petit. 2002. J. Clin. Microbiol. 40:2079-2083.
- Johnson, S., S.P. Sambol, J.S. Brazier, M. Delmée, V. Avesani, M.M. Merrigan, and D.N. Gerding. 2003. J. Clin. Microbiol. 41:1543-1547.
- Kato, H., N. Kato, K. Watanabe, N. Iwai, H. Nakamura, T. Yamamoto, K. Suzuki, S. Kim, Y. Chong, and E.B. Wasito. 1998. J. Clin. Microbiol. 36:2178-2182.
- Johnson, S., S.A. Kent, K.J. O'Leary, M.M. Merrigan, S.P. Sambol, L.R. Peterson, and D.N. Gerding. 2001. Ann. Intern. Med. 135:434-438.
- Wilkins, T.D. and D.M. Lyerly. 2003. J. Clin. Microbiol. 41:531-534.
- Fedorka, D.P. 2002. Clin. Microbiol. Newsletter. 24:76-79.
- Richmond, J.Y. and R.W. McKinney. 1999. Biosafety in Microbiological and Biomedical Laboratories. 4th ed. U.S. Government Printing Office. Washington, D.C.
- Lyerly, D.M., H.C. Krivan, and T.D. Wilkins. 1988. Clin. Microbiol. Rev. 1:1-18.
- Mandell, G.L., J.E. Bennett, and R. Dolin. 2000. Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases. 5th ed. Churchill Livingstone. New York, NY.

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Xpect® es una marca comercial de Remel Inc.